SQÛVA

SQ UV Activator

ENHANCE YOUR IMPLANT SURFACE

10 SECONDS IS ALL IT TAKES!



The Benefits of UV Irradiation!

DENTIS

6 Centerpointe Dr. Suite 600, La Palma CA 90623
Tel. 323-677-4363 | Fax. 323-677-4366 | info@dentisusa.com | www.dentisusa.com

Enhance Your Implant Surface
10 Seconds is All It Takes!

It directly makes the surface of the implant super hydrophilic after flashing the UV light.





Motion Sensor

Prevents cross contamination with contactless functionality

Actuator Control
- 1 second detection
  (LED Color: Green): Open/Close

UV processing Control
- 2 seconds detection
  (LED Color: Green ▶ Blue): Start

Voice Notifications

- Supports voice guidance during operation and when the LCD screen changes
- UV Processing Start, Complete

5" LCD Touch Panel

It's intuitive and can be simply controlled.

Setting Mode
- Filter Replacement, Sound
- Motion Sensor On/Off
- Language Choice (KOR/ENG)
- Sound Volume

User Convenience

O-Zone Sensor Installed
- Manage the user safety and prevents the accident through ozone leak detection

RFID Card Installed
- During product operation, its possible to insert and check the number of usages remaining and recharge

LED Effect
- Checking the operation status through lighting when the product is running

Benefits of SQUVA!

- ✓ Stable implant surgery by strengthening bone adhesion
- ✓ More satisfied patients with reduced surgery and healing time
- ✓ Gain trust with new patients with a higher standard of treatment

Why SQUVA?

- ✓ Can be placed immediately after extraction (When the bone defect is large)
- ✓ Can be used for elderly patients (When bone tissue and gum are in poor condition)
- ✓ Can be used on patients with chronic illness (Diabetes/high blood pressure)